

MyGo Five FAQs

Here are some frequently asked questions on the MyGo Five Kits.

FAQ	Answer
Do I need to freeze my MyGo Five on arrival?	No, it should be stored at room temperature until made up. Once hydrated, the kit should be stored at 4°C.
How long will my MyGo Five kit last for?	Your kit has a shelf life of 12 months from the date of shipment until the lyophilised reagents are hydrated. Once hydrated the reagents can be stored at 4°C for up to 12 hours. Hydrated reagents must be used within 12 hours of hydration.
How many runs will my MyGo Five give me?	You will get 48 reactions with your MyGo Five kit. This will be sufficient to perform 2 MyGo Five experiments on your MyGo Pro.
What disposables should I use with the kit?	On the MyGo Pro you should use MyGo 8 well strips and caps. On the MyGo Mini you should use MyGo single tubes.
What run parameters should I use with the kit?	You can download the experimental template from mygopcr.com . All the run parameters you will need will already be set up in this template.
What results should I expect?	You should expect 8 amplification curves per fluorescent dye. 4 of these are singleplex reactions and the other 4 will be multiplex reactions (blue).
Can I use the MyGo Five on other platforms?	No, the MyGo Five is only compatible on the MyGo Pro.